Safety Instructions / 安全提示 / Инструкции за безопасност
Bezpečnostní pokyny / Sicherheitshinweise / Sikkerhedshenvisninger
Ohutusjuhised / Indicaciones de seguridad / Consignes de sécurité
Υποδείξεις για την ασφάλεια / Sigurnosna napomena
Biztonsági utasítások / Treoracha sábháilteachta / Avvertenze di sicurezza
Saugumo nuorodos / Drošības norādījumi / Noti dwar is-sigurtà
Veiligheidswaarschuwingen / Wskazówki bezpieczeństwa
Instruções de segurança / Reguli de siguranță / Säkerhetsanvisningar
Turvallisuushjeet / Bezpečnostné pokyny / Varnostna opozorila

MK (E6) • MKF (E6)

www.binder-world.com/manuals

Download the full operating manual.

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Safety instructions

MK (E6)

Alternating climate chambers with program control

MKF (E6)

Alternating climate chambers with program control

Important note

The present document is only valid in conjunction with the operating manual that correspond to your device

- MK (E6) art. no. 7001-0714
- MKF (E6) art. no. 7001-0715

The initial commissioning and usage of the aforementioned devices (during the entire product life cycle) is only permitted with the operating manual. The operating manual must be read completely before usage.

The present document is not a replacement for the operating manual, as it only contains selected information.

You can download the appropriate operating instructions for your device using this QR code or at www.binder-world.com/manuals. If this is not possible, please contact BINDER GmbH:

Address	Post office box 102, 78502 Tuttlingen, Germany
Tel.	+49 7462 2005 0
Internet	www.binder-world.com
Service Hotline	+49 7462 2005 555
Service Fax	+49 7462 2005 93 555

Safety instructions

The following represents the hazards for people. Safety information with lower risk levels, and important detailed information about safety, can be found in the operating manual. To avoid injuries and damage observe all safety instructions of the operating manual. It is very important that the operating manual is read completely before initial commissioning and usage.

Pictograms

 Risk of injury. Observe the safety instructions in the operating manual	 Electrical hazard	 Hot surface	 Explosive atmosphere
 Stability hazard	 Lifting hazard	 Scalding hazard	 Risk of chemical burns
 Harmful substances	 Biohazard	 CO ₂ suffocation and poisoning hazard	
 Mandatory regulation	 Read operating instructions	 Disconnect the power plug	 Lift with mechanical assistance
 Wear protective gloves	 Wear safety goggles		
 Do NOT touch	 Do NOT spray with water	 Do NOT climb	

	 DANGER
	Electrical hazard during live maintenance work. Deadly electric shock. Ø Do NOT remove the rear panel of the chamber. ➤ Disconnect the chamber before conducting maintenance work. Turn off the main power switch and pull the power plug. ➤ Make sure that general maintenance work will be conducted by licensed electricians or experts authorized by BINDER.

	 DANGER Explosion or implosion hazard through the introduction of unsuitable loading material, or through combustible dusts or explosive mixtures in the vicinity of the chamber or in the chamber. Serious injury or death from burns and / or explosion pressure. - Ø Do NOT operate the chamber in potentially explosive areas. - Ø Make sure that there are NO combustible dust or air-solvent mixture in the vicinity of the chamber. - Ø Do NOT introduce any substance combustible or explosive at working temperature into the chamber, in particular no energy sources such as batteries or lithium-ion batteries - Ø Make sure that there are NO combustible dust or air-solvent mixture in the inner chamber.
	 DANGER Electrical hazard by water entering the chamber. Deadly electric shock. - Ø Do NOT install the chamber in damp areas or in puddles. ➢ Set up the chamber in a splash-proof manner. ➢ Before cleaning, turn off the chamber at the main power switch and disconnect the power plug. Let the chamber cool down to ambient temperature. - Ø The chamber must NOT become wet during operation, cleaning, or maintenance. Do NOT spill water or cleaning agents over the inner and outer chamber surfaces. ➢ Completely dry the chamber before turning it on again.
	 DANGER Electrical hazard through overload of contacts on the zero-voltage relay alarm output. Deadly electric shock. Damage to the switching contacts and connection socket. - Ø Do NOT exceed the maximum switching load of 24V AC/DC – 2.5A. - Ø Do NOT connect any devices with a higher loading capacity.
	 DANGER Electrical hazard due to missing protective conductor connection or high leakage current in case of improper installation. Deadly electric shock. ➢ Make sure that the chamber's power plug and the power socket match and securely connect the electrical protective conductors of the chamber and the house installation. ➢ Only connect the plug to a socket if it has an intact earth connection.

	 DANGER
	Risk of electric shock due to inadequate electrical protection. Fatal electric shock. Fire hazard. ➤ Use type B residual current devices. ➤ Ensure adequate electrical protection is provided for the number of chambers to be operated.

 	 WARNING
	Danger of intoxication and infection through contamination of the chamber with toxic, infectious or radioactive substances. Poisoning. Damages to health. - Ø Do NOT introduce any substance which could lead to release of toxic gases ➤ Protect the interior of the chamber from contamination by toxic, infectious or radioactive substances. ➤ Take suitable protective measures when introducing and removing toxic, infectious or radioactive material. - Ø NEVER take a chamber contaminated with toxic substances or sources of infection for recycling according to Directive 2012/19/EU. ➤ Prior to disposal, remove all toxic substances and sources of infection from the chamber. ➤ A chamber from which all toxic substances or sources of infection cannot be safely removed must be considered as "special" waste according to national law. Dispose of it accordingly.

  	 CAUTION
	Risk of injury and damages by lifting heavy loads and by sliding or tilting of the chamber due to improper lifting or transport. Injuries. ➤ For transport, secure the chamber with transport straps. ➤ Keep the chamber in upright position. ➤ Lift the chamber using technical devices (fork lifter). Set the fork lifter only from the rear in the middle of the chamber. Make sure to place all the lateral supports of the chamber on the forks. - Ø Do NOT load the lower housing cover with heavy objects while the chamber door is open and do NOT climb on it.

	 CAUTION
	Danger of burning by touching hot chamber parts during or after operation. Burns. Do NOT touch the inner surfaces, the front panel around the inner chamber, the window, the access port area or the charging material during operation.

	 CAUTION Danger of chemical burns through contact with skin or ingestion of the neutral cleaning agent or through eye contact with the disinfectant spray. Skin and eye damage. ➤ Do not ingest the neutral cleaning agent. Keep it away from food and beverages. ➤ Wear protective gloves and goggles. ➤ Avoid skin contact with the neutral cleaning agent.
---	--

	 DANGER Danger of suffocation and poisoning by high concentration of CO₂ (> 4 Vol.-%). Death by suffocation. Ø Do NOT set up chambers in non-ventilated recesses. ➤ Ensure technical ventilation measures. ➤ Observe the relevant regulations for handling CO₂. ➤ Ensure sufficient room volume. The required room volume must be calculated in accordance with EN 378-1:2016+A1:2020.
---	--

For MKF:

	 WARNING Risk of overheating or fire and risk of damage if the chamber continues to be operated with the alarm message “Humidity system”. Injuries and damage to the chamber and its surroundings. ➤ Do NOT continue to operate the chamber if the alarm message “Humidity system” appears. ➤ Do NOT acknowledge the alarm message “Humidity system”. ➤ If the alarm message “Humidity system” appears, turn off the chamber and contact BINDER service.
---	---

Intended use

The chambers are suitable for conditioning of harmless materials. They are suitable for temperature drying and heat treatment of solid or pulverized charging material, as well as bulk material, using the supply of heat. The chambers are specially designed for solving all the problems which occur during material and ageing tests.

Any solvent must not be explosive and flammable. A mixture of any component of the charging material with air must NOT be explosive. The operating temperature must lie below the flash point or below the sublimation point of the charging material. Any component of the charging material must NOT be able to release toxic gases. Other applications are NOT approved.

The chambers are designed for setting up inside a building (indoor use). The requirements described in the Operating Manual for installation site and ambient conditions must be met

The chambers are not classified as medical devices as defined by Regulation (EU) No 2017/745.

The loading material shall not contain any corrosive ingredients that may damage the machine components made of stainless steel, aluminum, and copper. Such ingredients include in particular acids and halides. Contamination of the chamber by toxic, infectious or radioactive substances must be prevented. The chambers do not dispose of any measures of explosion protection.

Details can be found in the operating manual in the chapter "Intended use". Any use of the chambers that does not comply with the requirements specified in the Operating Manual shall be considered improper use. Make sure that you have the Operating Manual that match your chamber, including any accessories.

You may not make any modifications to the chamber or device yourself, as these may limit proper usage.

To operate the chamber, use only original BINDER accessories or accessories / components from third-party suppliers authorized by BINDER. The user is responsible for any risk arising from using unauthorized accessories

BINDER GmbH is responsible for the safety features of the chamber only, provided skilled electricians or qualified personnel authorized by BINDER perform all maintenance and repair, and if components relating to chamber safety are replaced in the event of failure with original spare parts. The user is responsible for any risks arising from using unauthorized accessories/components.

For MKF:


Types of suitable water quality for freshwater supply:

- Deionized water from a water treatment installation already existing at the customer's site. Conductivity from 1 $\mu\text{S}/\text{cm}$ up to a maximum of 20 $\mu\text{S}/\text{cm}$. (Water, which is in equilibrium with the CO_2 in the air, and has a conductivity below 1 $\mu\text{S}/\text{cm}$ (ultrapure water), may cause acid corrosion due to its low pH)
- Water treated by the optional water treatment system BINDER Pure Aqua Service (disposable system).


Carbon dioxide (CO_2) is used as refrigerant. Vent out any CO_2 gas that may escape via good room ventilation or a suitable connection to an exhaust system. We recommend the installation of a CO_2 warning system. Ensure sufficient room volume. The calculation of the required room volume must be carried out according to EN 378-1:2016+A1:2020.

Chamber overview

Front


MK/MKF 56

MK/MKF 240

Operating panel on the left chamber side

- A Triangular instrument panel (MK/MKF 56, 115)
- B Triangular signal panel (MK/MKF 240, 720, 1020)
- C Instrument panel (MK/MKF 240, 720, 1020)
- D Inspection window
- E Door handle
- F Chamber door
- G MK: Refrigerating machine, maintenance access flaps
MKF: Refrigeration / humidity module
- H Main power switch
- I MKF: Access to fill the water can
- J Temperature safety device class 2 for over and under temperature (option): Entry displays for upper (J 1) and lower (J 2) temperature limit
- K Ethernet interface
- L RS485 interface (option)
- M 2 zero-voltage switching outputs via special controller functions (option for MK, regular equipment for MKF)
- N 2 zero-voltage switching outputs via special controller functions (option for MK, regular equipment for MKF)
- O Analog output for temperature (option)
- P Analog output for humidity (option for MKF)

Electrical connection data and installation

- ▶ The chambers come with a fixed power connection cable of at least 1800 mm / 70.87 in in length.
- ▶ Pollution degree acc. to IEC 61010-1: 2. Over-voltage category acc. to IEC 61010-1: II.
Housing protection acc. to EN 60529: IP 20. Internal over-current release category B 3 poles
- ▶ Nominal voltage and nominal power, 50 Hz units

	MK 240 9x20-0513	MK 400 9x20-0518	MK 720 9x20-0522	MK 1020 9x20-0527
Nominal voltage V (+/-10%) at 50 Hz	400	400	400	400
Current type	3N~	3N~	3N~	3N~
Nominal voltage kW	5.20	6.80	8.70	9.60
Nominal power A	8.90	11.90	14.00	22.30
Over-current release A	3x 16	3x 16	3x 16	3x 32
	MKF 240 9x20-0514	MKF 400 9x20-0519	MKF 720 9x20-0523	MKF 1020 9x20-0528
Nominal voltage V (+/-10%) at 50 Hz	400	400	400	400
Current type	3N~	3N~	3N~	3N~
Nominal voltage kW	6.80	7.90	11.00	11.50
Nominal power A	10.80	11.80	17.80	24.70
Over-current release A	3x 16	3x 16	3x 25	3x 32

- ▶ Nominal voltage and nominal power, 60 Hz units

	MK 240 9x20-0515	MK 400 9x20-0520	MK 720 9x20-0524	MK 1020 9x20-0529
Nominal voltage V (+/-10%) at 60 Hz	480	480	480	480
Current type	3N~	3N~	3N~	3N~
Nominal voltage kW	5.20	6.80	8.70	9.60
Nominal power A	7.50	10.00	13.50	17.60
Over-current release A	3x 16	3x 16	3x 16	3x 25
	MKF 240 9x20-0516	MKF 400 9x20-0521	MKF 720 9x20-0525	MKF 1020 9x20-0530
Nominal voltage V (+/-10%) at 60 Hz	480	480	480	480
Current type	3N~	3N~	3N~	3N~
Nominal voltage kW	6.80	7.90	11.00	11.50
Nominal power A	8.50	10.5	13.6	18.5
Over-current release A	3x 16	3x 16	3x 20	3x 25

- ▶ The domestic socket must also provide a protective conductor. Make sure that the connection of the protective conductor of the domestic installations to the chamber's protective conductor meets the latest technology. The protective conductors of the socket and plug must be compatible!
- ▶ When connecting, please observe the regulations specified by the local electricity supply company as well as the local or national electrical regulations.
- ▶ Observe a sufficient current protection according to the number of devices that you want to operate. Use residual current circuit breakers type B for this purpose.
- ▶ Prior to connection and start-up, check the power supply voltage. Compare the values to the specified data located on the chamber's type plate.
- ▶ To completely separate the chamber from the power supply, you must disconnect the power plug. Install the chamber in a way that the power plug is easily accessible and can be easily pulled in case of danger.

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