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**GB / US VP Pump Modules
for Vacuum Drying Ovens VD and VDL
Mounting Instructions**

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The products described above are in conformity with the following EU Directives:

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The products described above, corresponding to this, bear the CE-mark.

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According to the BINDER internal production control, the products described above only contain electrical components that meet the relevant standards or do not fall within their scope of application. The technical documents from the component manufacturers serve to prove conformity.

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**UK
CA**

UKCA Declaration of Conformity

Name and address of manufacturer	BINDER GmbH Im Mittleren Ösch 5, 78532 Tuttlingen, Germany
Name and address of UK Authorised Representative	Comply Express Ltd Unit C2, Coalport House, Stafford Park 1, Telford TF3 3BD
Object of the Declaration	Control cabinets for accommodating a vacuum pump for vacuum drying ovens VD
Type Designation	Pump modules for VD 23, 56, 115
BINDER Art. No.	8012-2718, 8012-2721, 8012-2719, 8012-2722, 8012-2720, 8012-2723

The Objects of the Declaration described above are in conformity with the relevant UK Regulations and UK Guidelines:

- **Electrical Equipment (Safety) Regulations 2016**
Statutory Instruments 2016 No. 1101 – Consumer Protection Health and safety
- **Electromagnetic Compatibility Regulations 2016**
Statutory Instruments 2016 No. 1091 – Electromagnetic Compatibility
- **The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012**
Statutory Instruments 2012 No. 3032 – Environmental Protection

References of standards and/or technical specifications applied for this Declaration of Conformity, or parts thereof:

S.I. 2016 No. 1101:	According to the BINDER internal production control, the products described above only contain electrical components that meet the relevant standards or do not fall within their scope of application. The technical documents from the component manufacturers serve to prove conformity.
S.I. 2016 No. 1091:	EN IEC 61326-1:2021
S.I. 2012 No. 3032:	EN IEC 63000:2018

This Declaration is issued under the sole responsibility of the manufacturer.

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Mounting instructions for pump module

Contents

1. Safety	35
1.1. Personnel Qualification	35
1.2. Mounting Instructions	35
1.3. Legal considerations	35
1.3.1. Intellectual property	36
1.4. Structure of the safety instructions.....	36
1.4.1. Signal word panel	36
1.4.2. Safety alert symbol	37
1.4.3. Pictograms	37
1.4.4. Word message panel structure	37
1.5. Localization / position of safety and information labels on the chamber	38
1.6. Type plate of the VD pump modules.....	38
1.7. UKCA Label	39
2. Scope of delivery	40
2.1. Included connection parts	40
2.2. Accessories for VP1 / VP2	42
2.3. Accessories for VP5 EX.....	43
3. Mounting the vacuum drying oven on the pump module.....	44
3.1. Installation	44
3.2. Mounting the connecting plate	45
4. Installing the suction line	46
4.1. Fixing the hose at the housing rear.....	46
4.2. Hose connection at the vacuum pump VP1 / VP2	47
4.3. Hose connection at the vacuum pump VP5 EX.....	48
5. VDL: Inserting / removing the condensate collecting tray	49
6. Achieving equipotential bonding	50
6.1. Equipotential bonding between the vacuum drying oven and pump module	50
6.2. VDL: Equipotential bonding between the condensate collecting tray and the pump module	50
6.3. Equipotential bonding between the vacuum pump VP5 EX and the pump module	52
7. Placing the VP1 / VP2 vacuum pump in the VD pump module	54
7.1. Positioning the pump.....	54
7.2. Ventilation concept.....	55

8. **Placing the VP5 EX vacuum pump in the VDL pump module57**

9. **VDL: Equipotential bonding between the pump module and the conductive surface of the installation site58**

10. **Connection of a suction system at the pump module59**

10.1. Suction system at the VD pump module59

10.2. Suction system at the VDL pump module59

1. Safety

1.1. Personnel Qualification

The device must only be installed, tested, and started up by personnel qualified for assembly, startup, and operation of the device. Qualified personnel are persons whose professional education, knowledge, experience and knowledge of relevant standards allow them to assess, carry out, and identify any potential hazards in the work assigned to them. They must have been trained and instructed, and be authorized, to work on the device.

The device should only be operated by laboratory personnel especially trained for this purpose and familiar with all precautionary measures required for working in a laboratory. Observe the national regulations on minimum age of laboratory personnel.

1.2. Mounting Instructions

These Mounting Instructions are part of the scope of delivery. Always keep them handy for reference in the vicinity of the device. If selling the device, hand over the operating manual to the purchaser.

To avoid injuries and damage, observe the safety instructions given with the instructions.

	 CAUTION Dangers due to failure to observe the instructions and safety precautions. Serious injuries and damage to the device. ➤ Observe the safety instructions in these Mounting Instructions. ➤ Follow the operating procedures in these Mounting Instructions. ➤ Carefully read the complete Mounting Instructions of the device prior to installing and using it. ➤ Keep the Mounting Instructions for future reference.
--	--

	Make sure that all persons who use the device and its associated work equipment have read and understood the Mounting Instructions.
---	---

The Mounting Instructions are supplemented and updated as needed. Always use the most recent version of the manual. When in doubt, call the BINDER Service Hotline for information on the up-to-dateness and validity of these Mounting Instructions.

1.3. Legal considerations

These Mounting Instructions are for informational purposes only. They contain information for correct and safe installing, start-up, operation and decommissioning of the product. Note: the contents and the product described are subject to change without notice.

Understanding and observing the instructions in these Mounting Instructions are prerequisites for hazard-free use and safety during operation and maintenance. Images are to provide basic understanding. They may deviate from the actual version of the device. The actual scope of delivery can, due to optional or special design, or due to recent technical changes, deviate from the information and illustrations in these instructions. In no event shall BINDER be held liable for any damages, direct or incidental arising out of or related to the use of these instructions.

These Mounting Instructions cannot cover all conceivable applications. If you would like additional information, or if special problems arise that are not sufficiently addressed in this manual, please ask your dealer or contact us directly, e.g. by phone at the number located on page one of this manual.

Furthermore, we emphasize that the contents of these Mounting Instructions are not part of an earlier or existing agreement, description, or legal relationship, nor do they modify such a relationship. All obligations on the part of BINDER derive from the respective purchase contract, which also contains the entire and exclusively valid statement of warranty administration and the general terms and conditions, as well as the legal regulations valid at the time the contract is concluded. The statements in this manual neither augment nor restrict the contractual warranty provisions.

1.3.1. Intellectual property

These Mounting Instructions are protected by copyright. Any unauthorized copying or disclosure to third parties is strictly prohibited. We reserve the right to take legal action and, if necessary, to assert claims for damages in the event of infringement.

Trademark Information: All BINDER trademarks relating to products or service, as well as trade names, logos and product names used on the website, products and documents of BINDER company are trademarks or registered trademarks of BINDER company (including BINDER GmbH, BINDER Inc.) in the U.S. and other countries and communities of states. This includes word marks, position marks, word/figurative marks, design configurations, figurative marks, and design patents.

Patent Information: BINDER products, categories of products, and accessories may be covered by one or more patents and/or utility models in the U.S. and other countries and communities of states. This information is provided to satisfy the virtual patent marking provisions of various jurisdictions, in particular it is intended to serve as notice under 35 U.S.C. § 287(a). Products and services listed on the BINDER website may be sold individually or as part of a combination product. Additional patent applications may also be pending in the U.S. and other countries and communities of states.

Please visit www.binder-world.com for more information.

1.4. Structure of the safety instructions

In these Mounting Instructions, the following harmonized denominations and symbols indicate dangerous situations following the harmonization of ISO 3864-2 and ANSI Z535.6.

1.4.1. Signal word panel

Depending on the seriousness and probability of the consequences, dangers are identified with a signal word, the corresponding safety color, and if appropriate, the safety alert symbol.

 WARNING
Indicates an imminently hazardous situation that, if not avoided, will result in death or serious (irreversible) injury.

 CAUTION
Indicates a potentially hazardous situation which, if not avoided, may result in moderate or minor (reversible) injury.

NOTICE

Indicates a potentially hazardous situation which, if not avoided, may result in damage to the product and/or its functions or of a property in its proximity.

1.4.2. Safety alert symbol

Use of the safety alert symbol indicates **risk of injury**.

Observe all measures that are marked with the safety alert symbol in order to avoid death or injury.

1.4.3. Pictograms

Warning signs			
			
Electrical hazard	Danger of frost	Harmful substances	Pressurized device
			
Explosive atmosphere	Stability hazard	Lifting hazard	
Mandatory action signs			
			
Read instructions	Environment protection	Lift with several persons	

1.4.4. Word message panel structure

Type / cause of hazard.

Possible consequences.

- ⊘ Instruction how to avoid the hazard: prohibition
- Instruction how to avoid the hazard: mandatory action

Observe all other notes and information not necessarily emphasized in the same way, in order to avoid disruptions which could result in direct or indirect injury or property damage.

1.5. Localization / position of safety and information labels on the chamber

The following labels are located on the device:

Pictograms (Warning signs)	
	Hot surface (on the outer door of the pump modules)
	Risk of injury. Observe the safety instructions in the operating manual. (next to the socket of the VD pump modules)
	Keep safety labels complete and legible.

Replace safety labels that are no longer legible. Contact BINDER Service for these replacements.

1.6. Type plate of the VD pump modules

The type plate sticks to the left side of the pump module, bottom right-hand.

Nominal temp.	1,80 kW / 8,0 A		
	200-230 V / 50 Hz		
IP protection	20		
Safety device	200-230 V / 60 Hz		
Class	1 N ~		
Art. No.	8012-2718		
Project No.			
Built	2025	Pump Module VD	
		BINDER GmbH Im Mittleren Ösch 5 78532 Tuttlingen / Germany www.binder-world.com	
		VD-M 23	Serial No. 00000000000000
			Made in Germany

Type plate (example pump module for VD 23, 230 V)

Indications of the type plate

Indication		Information
BINDER		Manufacturer: BINDER GmbH
VD-M 23		Model designation
Serial No.	00000000000000	Serial no.
Nominal temperature		---
Enclosure protection	IP 20	IP type of protection acc. to standard EN 60529
Temp. safety device	DIN 12880	---
Class	---	---
Art. No.	8012-2579	Art. no.
Project No.		

Indication	Information
0,02 kW	Nominal power
200-230 V, 50 Hz 200-230 V, 60 Hz	Nominal voltage range +/- 10% at the indicated power frequency
1 N ~	Current type
8,0 A	Nominal power

Symbols on the type plate

Symbol	Information
	CE conformity marking
	Electrical and electronic equipment manufactured / placed on the market in the EU after 13 August 2005 and be disposed of in separate collection according to Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).

1.7. UKCA Label

The sticker with UKCA Authorised Representative details sticks next to the type plate to the left side of the pump module, bottom right-hand.

UKCA Label

Symbol on the sticker

Symbol	Information
	UKCA conformity marking

2. Scope of delivery

2.1. Included connection parts

To place the vacuum drying oven on top of the pump module

Connection parts, included with the pump module

3 x Brackets with cable ties


Connecting plate


To establish equipotential bonding

Grounding cables, included with the pump module

VP5 EX: Grounding cables for connection of the condensate collection tray and of the pump (pre-installed)


VP5 EX: Grounding cable for connection to conductive surface (pre-installed)


VP1 / VP2:


Ventilation sleeve for pump motor Art. no. 6016-0070, included with the pump

Upper silicone sleeve, for installation see chap. 7.2



VP5 EX:

Condensate collection tray Art. no. 4022-0495, included with the pump module	
Condensate collection tray (pre-installed)	


To connect the VP1 / VP2 vacuum pumps to VD series vacuum drying ovens

Connection kit VD VP1-VP2 Art. no. 8012-2613	
Small flange with hose nozzle, small flange DN 16 / DN 10, PP	
Aluminum straining ring, small flange DN 10/16 for connection on the VD	
Outer centering ring, small flange DN 10/16, with sealing ring NBR	
Elastomer vacuum hose DN 10, length 1900 mm	

To connect the VP5 EX vacuum pump to VDL series vacuum drying ovens

Connection kit VDL VP5 stainless steel Art. no. 8012-2724

Aluminum straining ring, small flange DN 10/16 for connection on the VDL	
Outer centering ring, small flange DN 10/16, with sealing ring NBR	
Metal hose, stainless steel 1.4404 / AISI 316, small flange DN 16, length 2000 mm	

2.2. Accessories for VP1 / VP2

Emission condenser VP1-VP2 Art. no. 8012-2615

Mounting kit with exhaust-side exhaust condenser. Contains all necessary small parts.	
---	---

Please note the additional information and installation instructions in the supplementary sheet for the emission condenser VP1-VP2, which is included with the article.

Membrane wrench SW 66 Art. no. 8500-0542

Special tool for replacing the membranes	
--	---

Service kit VP1 Art. no. 8500-0536

Wear parts set, consisting of membranes, valves, and seals (O-rings)	
--	--

Service kit VP2 Art. no. 8500-0537

Wear parts set, consisting of membranes, valves, and seals (O-rings)

2.3. Accessories for VP5 EX

To connect the VP5 EX vacuum pump to VDL series vacuum drying ovens

Connection kit VDL VP5 PTFE Art. no. 8012-2614

2 x Clamping ring, aluminum, small flange, KF DN 10/16 for connection on VDL

2 x External centering ring, small flange DN 10/16, with sealing ring NBR

PTFE antistatic tubing, small flange DN 16/1900 mm

Please note the additional information and installation instructions in the supplementary sheet for the VDL VP5 PTFE connection kit, which is included with the connection kit.

Service kit VP5 EX Art. no. 8500-0538

Wear parts set, consisting of membranes, valves, and seals (O-rings)

Membrane wrench SW 66 Art. no. 8500-0542

Special tool for replacing the membranes

Sealing ring Art. no. 8500-0540

O-ring 28 x 2.5
on the spherical ground joint of the round-bottom flask


3. Mounting the vacuum drying oven on the pump module

3.1. Installation

Place the vacuum drying oven on top of the pump module.

	 CAUTION Sliding, falling down, or tilting of the unit. Damage to the unit. Risk of injury by lifting heavy loads. - Ø Do NOT lift the unit using the door handle, the door or at the lower housing. ➤ Lift chambers size 23 and 56 at its four lower corners with the aid of four people. ➤ Lift the chambers size 115 with the aid of six people.
---	---

On the pump module top surface there are two protruding hexagon head screws. The feet of the vacuum drying oven provide a central recess for these screws. Place the corner feet exactly on these screws. Thus, the chamber is protected against lateral sliding.


3.2. Mounting the connecting plate

An additional consolidation is obtained by mounting the delivered connecting plate. Unscrew the according Phillips screws at the vacuum drying oven and at the pump module and use them to fasten the connecting plate.

Fixation of the connecting plate:

Sizes 23 and 56: Fixation with 6 screws

Size 115: Fixation with 8 screws

Mounting this connecting plate is part of the equipotential bonding concept.


4. Installing the suction line

4.1. Fixing the hose at the housing rear

Before attaching the hose to the rear wall of the vacuum drying oven, it must be connected to the hose connection of the vacuum drying oven.

Insert the supplied cable connector into the 3 provided holes at the housing rear of the vacuum drying oven and fasten the hose with the cable ties.

Install the tube downgrade toward the pump and conduct it from behind into the pump module.


Vacuum drying oven with mounted pump module and connected hose

4.2. Hose connection at the vacuum pump VP1 / VP2

With the accessory “Pump module with pump”, the pump module comes with the VD VP1-VP2 connection kit (Art. No. 8012-2613), which includes besides an appropriate vacuum hose all necessary connection parts to connect the hose and the vacuum pump.

Connect the vacuum hose to the suction connection of the pump. To do this, plug the vacuum hose DN 10 onto the hose shaft at the pump inlet.

4.3. Hose connection at the vacuum pump VP5 EX

With accessory “Pump module with pump”, the pump module comes with the connection kit VDL VP5 stainless steel (Art. no. 8012-2724). The connection kit VDL VP5 PTFE (Art. no) is optionally available. Each connection kit includes an appropriate vacuum hose and all necessary connection parts to connect the hose and the vacuum pump.


Place the vacuum pump in the pump module (chap. 8).

Establish the equipotential bonding: To create a conductive connection between the vacuum pump and the pump module, a grounding cable is supplied that is already connected to the pump module. Attach it to the pump using the material supplied (grounding kit). (chap. 6.3).

Connect the connection kit to the pump as described in the manual of the VDL VP5 PTFE connection kit.

Please note the information and installation instructions in the supplementary sheet to the PTFE connection kit that is included with the article.

5. VDL: Inserting / removing the condensate collecting tray


The condensate collection tray for the VDL pump module is delivered fully assembled, including grounding.

Safety concept: If solvent is accidentally spilled when removing the filled condensate flask from the pump, it will be collected in the tray. Additional safety measures include suction from the pump module (chap. 10) and the prescribed installation of the pump module on a conductive surface with a conductive connection between the pump module, pump, and collection tray (chap. 9).

The condensate collection tray can be removed for **cleaning and emptying**; it is simply plugged in.

Removal: Lift the condensate collection tray at the front of the frame and then carefully remove it. Be careful not to spill any condensate. Make sure you wear conductive clothing/personal protective equipment.

Positioning: Insert the condensate collection tray diagonally from top to back into the pump module and then insert the guide latch at the front center into the slot in the pump module (not shown).

Procedure:

- Place the condensate collection tray in the pump module as described.
- Then establish the connection to the equipotential bonding system as described in chap. 6.2.

6. Achieving equipotential bonding

Equipotential bonding must always be made via external grounding connections, so that no potential can be introduced in the event of a short circuit. Following the installation of the vacuum drying oven, the pump module, and the pump and implementing all the measures described below to achieve equipotential bonding, we recommend performing a protective conductor measurement before commissioning.


6.1. Equipotential bonding between the vacuum drying oven and pump module

The conductive connection between the vacuum drying oven and the pump module is established by mounting the connecting plate (chap. 3.2). Therefore, in addition to the stabilizing function, this plate also fulfills a protective function by ensuring equipotential bonding


6.2. VDL: Equipotential bonding between the condensate collecting tray and the pump module

Delivery condition: The condensate collection tray is already installed and the grounding cable is plugged in at the factory. Check that the cable lug is securely in place.

To empty the tray, the grounding cable must be removed from the cable lug and reconnected after the tray has been reinserted.


View into the opened pump module with the two grounding cables for the condensate collecting tray and vacuum pump




Position of the stud bolt back right at the bottom of the pump module (example size 23) behind the condensate collection tray with connections of both grounding cables for the condensate collection tray (bottom) and vacuum pump (top)
(view from behind, rear wall of the pump module not shown)

Procedure:

- Insert the condensate collection tray into the pump module (chap. 5)
- To establish the conductive connection, a equipotential bonding cable (green-yellow, cross-section 6 mm²) with a cable lug at the free end is supplied. This cable is already attached to the stud bolt at the rear right of the pump module bottom.
- Plug the cable lug of the grounding cable onto the grounding lug of the condensate collection tray.


6.3. Equipotential bonding between the vacuum pump VP5 EX and the pump module

	View inside the open pump module showing with both grounding cables for the condensate collection tray and the vacuum pump
	Position of the stud bolt back right at the bottom of the pump module (example size 23) behind the condensate collection tray with connections of both grounding cables for the condensate collection tray (bottom) and vacuum pump (top) (view from behind, rear wall of the pump module not shown)
	stud bolt back right at the bottom of the pump module with connections of both grounding cables for the condensate collection tray (bottom) and vacuum pump (top)


Procedure:

- To establish a conductive connection between the vacuum pump and the pump module a grounding cable (green-yellow, diameter 6 mm²) with ring cable lug is supplied. This cable is already fixed at the stud bolt back right at the bottom of the pump module. The other fastening material is already pre-installed.
- Use the antistatic connection point on the pump housing to attach the equipotential bonding cable to the pump module
- No existing ground connection may be loosened or removed!

Afterwards, place the pump correctly in the pump module (chap. 8).

Connection point for the antistatic connection


View of the pump from behind after making the equipotential bonding connections and placing the pump in the pump module

7. Placing the VP1 / VP2 vacuum pump in the VD pump module

7.1. Positioning the pump



To ensure cooling of the vacuum pump the pump must be positioned precisely. The pump feet are used for centering. Additionally, rubber buffers are attached to the floor of the pump module.

Holes in the floor of the pump module for the pump feet, and rubber buffers


Procedure:

- Install the vacuum pump in a way that its four feet sit in the provided holes of the pump module bottom.


Placement of a vacuum pump foot in one of the holes in the floor of the pump module

After the pump is placed, the weight of the pump compresses the buffers.

7.2. Ventilation concept


The motor is ventilated through a hole in the bottom of the pump module. To do this, the pump must be positioned precisely.


Representation of the air flow

Two silicone ventilation sleeves are used to ensure the pump's motor cooling.

The lower silicone sleeve is already installed; the upper silicone sleeve, which is included with the pump module, must be installed:



lower silicone sleeve, pre-installed

upper silicone sleeve, after installation

Procedure:

- Clamp the upper ventilation sleeve above the base plate of the pump as shown:

Vacuum pump with both silicone sleeves after insertion


This means that the pump sits with its fan grille exactly above the hole and sucks in the air from outside, not from inside the pump module.

8. Placing the VP5 EX vacuum pump in the VDL pump module

Before installing the vacuum pump, equipotential bonding with the pump module must be established (chap. 6.3).

The pump feet are used for correct positioning.

Square recesses in the condensate collection tray form recesses for the pump feet.

Placing the VP5 EX vacuum pump
in the condensate collection tray
(example volume 23)

9. VDL: Equipotential bonding between the pump module and the conductive surface of the installation site

To establish a conductive connection between the pump module and the conductive surface of the installation site a grounding cable (green-yellow, diameter 6 mm²) with ring cable lugs is supplied with the pump module.


	
Fixation screw fitting on the rear of the pump module	Grounding cable, mounted (state of delivery)

Procedure:

- If necessary: Attach the grounding cable to the housing of the pump module. For this purpose the pump module, like the VDL, provides a fixation screw fitting for grounding on the bottom left of the rear panel. Unscrew the nut and fix the ring cable lug of the grounding cable on stud bolt of the screw fitting. Secure the connection with a toothed disk and screw on the nut again.
- Connect the other end of the grounding cable with the conductive surface of the installation site.

10. Connection of a suction system at the pump module

On the rear wall of the pump module there is a connection of 100 mm diameter, approx. 40 mm in length.

10.1. Suction system at the VD pump module

You can optionally connect a suction system to the provided exhaust air connection on the rear panel of the VD pump module. From a functional perspective, this is not necessary.

When using an optional suction system, we recommend a volume flow of 165 m³/h to ensure proper pump function.

10.2. Suction system at the VDL pump module

Connect a suction system to the provided exhaust air connection on the rear panel of the VDL pump module with a specific minimum flow rate. This results from the interaction of the set-point temperature and the ambient temperature and is determined by the operator.

In addition, the suction system serves to prevent the formation of an explosive atmosphere inside the pump module in case any solvents are accidentally spilled when removing the filled condensate collector flask

In addition, the entire installation area should be provided with a technical ventilation (suction), which serves to prevent formation of an explosive atmosphere during loading of the oven as well as during removal of the charging material after incomplete drying.

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