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Installation Manual for the Ethernet/RS422-RS485 Converter

Content

1.	Information notice regarding the manual.....	12
2.	Intended use.....	12
3.	Unit description.....	13
4.	Scope of delivery of Ethernet/RS422-RS485 interface converter kit	14
5.	Possibilities of networking.....	15
5.1.	Each temperature chamber via its own Ethernet/RS422-RS485 converter	15
5.2.	All temperature chambers via the same Ethernet/RS422-RS485 converter	16
5.3.	Bus termination / RS485 operation	16
6.	Connection of the Ethernet/RS422-RS485 converter to the Ethernet.....	17
7.	Logical connection of the MAC address to the IP address	18
8.	Driver installation for the Ethernet/RS422-RS485 converter for Windows NT4, 2000, XP 20	
9.	Service	21

1. Information notice regarding the manual

This installation manual contains a short description how to connect and install the Ethernet/RS422-RS485 converter and the drivers needed. It is expected that the converter has been purchased at the BINDER GmbH and thus the correct internal settings have already been carried out. It is also assumed that only devices and cables of the BINDER GmbH offered specially for this purpose will be used. To allow also different kinds of use offered by the Ethernet/RS422-RS485 converter, we include next to this manual also the original operation manual of the producer W & T. We recommend that the described procedures will be carried out by the customer's network administrator.

2. Intended use

The Ethernet/RS422-RS485 converter is suitable to network temperature chambers with RS 422 interface or RS485 interface, measurement devices, and the Alarm Box AB01 of the BINDER GmbH via an existing Ethernet network.

Thus networking is possible without additional laying of cables inside existing company-owned Ethernet networks. At any network computer the communication software APT-COM™ of the BINDER GmbH or an existing data acquisition system can be used to control the temperature chambers connected to the Ethernet and to receive and save the chamber measurement data.

3. Unit description

The Ethernet/RS422-RS485 converter is independent electrical equipment with external power supply via a 24V AC power supply unit. It can be loosely installed or mounted to clamping profile bars in distributor boxes. By internal DIP switches and plug bridges different operational modes can be set. The settings necessary for operating temperature chambers of the BINDER GmbH have already been pre-set by the BINDER GmbH. The connection to an Ethernet network (TCP/IP protocol) is established using a RJ45 socket and a PATCH cable. The connection to a temperature chamber of the BINDER GmbH or to a plug distributor for connecting several temperature chambers is established using a 9-poles Sub-D socket and a special cable. Both connection cables and the power supply unit are included with the scope of delivery.

The Ethernet/RS422-RS485 converter provides a fix MAC address to which an also fix IP address must be assigned in the Ethernet. A COM interface is assigned to the measurement and evaluation computer on which e.g., the communication software APT-COM™ of the BINDER GmbH is running, using a driver that will be installed on this computer. The Ethernet/RS422-RS485 converter is thus functionally transformed into a COM interface for this computer, even though it might be located at a far distinct place.

The connected temperature chamber(s) are addressed via this COM interface and the address that can be individually adjusted at each chamber controller.

Figure 1: Connected converter

4. Scope of delivery of Ethernet/RS422-RS485 interface converter kit

Complete kit:

Description	Art. no.
Set Schnittstellenwandler Ethernet/RS422-RS485	8012-0380

Components:

Description	Art. no.
Ethernet/RS422-RS485 converter with original description of producer W&T (only needed in case of problems during start-up)	5021-0016
Patch cable 3 m for connection to Ethernet	5023-0097
Connection cable 5 m for connecting the converter to an individual chamber or to a plug distributor RS422 from BINDER	5023-0095
Anleitung Anschluss Schnittstellenwandler Ethernet/RS422-RS485	7001-0068


Optional:

Description	Art. no.
Extension cable 50m (5023-0117) for connecting the connection cable 5m to an individual chamber or to a plug distributor RS422 from BINDER	5023-0117

5. Possibilities of networking

Depending on the local situation one of the following possibilities can be selected to connect the controlling and evaluation computer to the temperature chambers.


5.1. Each temperature chamber via its own Ethernet/RS422-RS485 converter

Notes: If the Ethernet/RS422-RS485 converter is already built-in in the temperature chambers (optional), only the PATCH cable will be plugged into the Ethernet interface (RJ45 socket) of the chamber. In this case the connection to the chamber controller is already established in a similar way as described here.

The optional extension cable 50m (5023-0117) can be connected behind the 5m cable (5023-0095).

5.2. All temperature chambers via the same Ethernet/RS422-RS485 converter

The optional extension cable 50m (5023-0117) can be connected behind the 5m cable (5023-0095).

5.3. Bus termination / RS485 operation

For the connection of RS485 networks, DIP switches 1 and 2 must be set to ON

To connect temperature chambers from BINDER, the recommended operating mode of the interface converter is termination enabled (BINDER factory setting).

Some controllers, like RD3, MB1, MB2, have an internal termination. When combining different controllers, communication errors can occur which are indicated by APT-COM™ as “timeout” or “CRC error”. In this case it can be necessary to disable termination at the interface converter.

To disable termination, set DIP switches 6 and 7 in the interface converter to Off. You can safely try out the correct setting. No damage of controller or interface converter will occur.


Figure 2: Termination setting via DIP switches 6 and 7

6. Connection of the Ethernet/RS422-RS485 converter to the Ethernet

Note: When connecting and installing the Ethernet/RS422-RS485 converter, it is required to consult the responsible network administrator. The BINDER service hotline is available for questions about the installation, provided that the Ethernet/RS422-RS485 converter has been purchased via the BINDER GmbH.

Connect the 24V AC power supply unit to the Ethernet/RS422-RS485 converter. (The plug-type connection is protected against reverse supply connection). Plug-in the power supply unit to a wall outlet appropriate to the connection data of the power supply unit's type plate. The green LED (POWER) at the Ethernet/RS422-RS485 converter must be illuminated. The red Error LED flashes.

Use the PATCH cable (3 m) included in the scope of delivery to connect the socket of the Ethernet/RS422-RS485 converter indicated as "NETWORK" to an RJ 45 outlet of the Ethernet network close to the temperature chamber(s). If the network connection is correct, the red LED must not flash.

7. Logical connection of the MAC address to the IP address

If the Ethernet/RS422-RS485 converter is connected to the Ethernet, the internally preset fix MAC address must be connected to a fix IP address of the Ethernet. The network administrator must assign a fix IP address to each Ethernet/RS422-RS485 converter. This IP address must not change in the Ethernet and must not be assigned more than once.

The connection of the MAC address to the IP address must be carried out individually for each Ethernet/RS422-RS485 converter. For this purpose the MS-DOS Command Prompt under Windows is started on any computer connected to the Ethernet.

For Windows 2000 and NT: Start → Programs → Accessories → Command Prompt

For Windows 95, 98, ME: Start → Programs → MS-DOS Command Prompt

The following command line establishes the logical connection of the Ethernet/RS422-RS485 converter's MAC address to the fix IP address assigned to it:

Q: \	>arp -s	192.168.10.38	00-c0-3d-01-35-61
	command	assigned IP address (example)	MAC address (example)

Notes:

Blanks are symbolized in the above example by underlined blanks.

The signs preceding the Command Prompt (Q:\>) vary in every network.

The preset fix MAC address of each Ethernet/RS422-RS485 converter is noted on the housing of the converter.



Figure 3: Entry of logical connection

If by this command this logical connection has been successfully transferred to the Ethernet/RS422-RS485 converter, the latter can be reached from now on in the Ethernet using its fix IP address.

Addressing the Ethernet/RS422-RS485 converter by its assigned IP address can be tested with the following command after the MS-DOS Command Prompt:

Q: \	>ping	192.168.10.38
	command	assigned IP address (example)

Correct addressing behavior is displayed as follows:


```

C:\> Command Prompt
Microsoft Windows 2000 [Version 5.00.2195]
(C) Copyright 1985-1999 Microsoft Corp.

D:\>ping 192.168.10.10

Pinging 192.168.10.10 with 32 bytes of data:

Reply from 192.168.10.10: bytes=32 time<10ms TTL=64
Reply from 192.168.10.10: bytes=32 time<10ms TTL=64
Reply from 192.168.10.10: bytes=32 time<10ms TTL=64
Reply from 192.168.10.10: bytes=32 time<10ms TTL=64

Ping statistics for 192.168.10.10:
    Packets: Sent = 4, Received = 4, Lost = 0 (0% loss),
    Approximate round trip times in milli-seconds:
        Minimum = 0ms, Maximum = 0ms, Average = 0ms
  
```

Figure 4: Positive result when checking the addressing behavior

Incorrect addressing behavior is displayed as follows:

```

C:\> Command Prompt
Microsoft Windows 2000 [Version 5.00.2195]
(C) Copyright 1985-1999 Microsoft Corp.

D:\>ping 192.168.10.38

Pinging 192.168.10.38 with 32 bytes of data:

Request timed out.
Request timed out.
Request timed out.
Request timed out.

Ping statistics for 192.168.10.38:
    Packets: Sent = 4, Received = 0, Lost = 4 (100% loss),
    Approximate round trip times in milli-seconds:
        Minimum = 0ms, Maximum = 0ms, Average = 0ms
  
```

Figure 5: Negative result when checking the addressing behavior

Some possible error causes:

The IP address has not been correctly or fix assigned in the network, or it has not been assigned to the network connection socket to which the Ethernet/RS422-RS485 converter is connected via the PATCH cable.

The power supply unit of the Ethernet/RS422-RS485 converter is not connected.

8. Driver installation for the Ethernet/RS422-RS485 converter for Windows NT4, 2000, XP

The last step of the Ethernet/RS422-RS485 converter installation is to establish a logical connection of the previously assigned fix IP address to a COM port of the controlling and evaluation computer on which APT-COM™ or other customer data acquisition software is running.

For this the driver delivered on CD is needed.

The driver is installed on the computer through an installation set-up. Start the file "Setup.exe" in the following directories, depending on the operating system.

For Windows 2000, X, and NT4: Directory "ComServ_Win2000_XP_NT4"

For Windows 95,98, and ME: Directory "ComServ_Win95_98_ME"

The set-up automatically leads through all installation steps.

The Setup installs the icon "COM Port Redirector" in the Windows Control Panel.

When clicking on this icon a dialogue window "W&T COM Port Redirector" opens for further driver settings.


Figure 6: Dialogue window for driver settings

By clicking on the button "COM Server" another dialogue window opens, in which the connection between a free computer COM port and the previously assigned fix IP address can be entered. The first free computer COM port is shown, so just enter the IP address. To do this, remove the check at "Set target by Registry ". If necessary also a different free computer COM port can be entered, for this press the button "Port Number". The other settings generally do not need to be altered.


Figure 7: Example of connection of an IP address to the COM-Port under Windows 2000

The Ethernet/RS422-RS485 converter is now completely installed. The temperature chambers are now registered in the communication software APT-COM™ with the COM number that has been assigned to the Ethernet / RS 422 converter they are (each or altogether, depending on the way of networking, see chap. 5.1 and 5.2) connected to. The chamber address required for the unambiguous chamber identification must be entered by hand individually and only once at the controller of each temperature chamber or measurement device. For a description see the respective operating manual of the chamber or device bearing an RS 422 interface or RS485 interface.

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

```
Command Prompt
Microsoft Windows 2000 [Version 5.00.2195]
(C) Copyright 1985-1999 Microsoft Corp.

D:\>ping 192.168.10.10

Pinging 192.168.10.10 with 32 bytes of data:

Reply from 192.168.10.10: bytes=32 time<10ms TTL=64
Reply from 192.168.10.10: bytes=32 time<10ms TTL=64
Reply from 192.168.10.10: bytes=32 time<10ms TTL=64
Reply from 192.168.10.10: bytes=32 time<10ms TTL=64

Ping statistics for 192.168.10.10:
    Packets: Sent = 4, Received = 4, Lost = 0 (0% loss),
    Approximate round trip times in milli-seconds:
        Minimum = 0ms, Maximum = 0ms, Average = 0ms
```

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```
Command Prompt
Microsoft Windows 2000 [Version 5.00.2195]
(C) Copyright 1985-1999 Microsoft Corp.

D:\>ping 192.168.10.38

Pinging 192.168.10.38 with 32 bytes of data:

Request timed out.
Request timed out.
Request timed out.
Request timed out.

Ping statistics for 192.168.10.38:
    Packets: Sent = 4, Received = 0, Lost = 4 (100% loss),
    Approximate round trip times in milli-seconds:
        Minimum = 0ms, Maximum = 0ms, Average = 0ms
```

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