BINDER

PURE AQUA SERVICE
(E2) Art.No.8012-0759

D	**Wasseraufbereitungssystem** **Montage- und Betriebsanleitung**
GB / US	**Water treatment system** **Mounting and Operating Instructions**
F	**Système de traitement des eaux** **Manuel de montage / Mode d'emploi**
ES	**Sistema de tratamiento de agua** **Manual de montaje y de funcionamiento**
RU	**Система очистки воды** **Инструкции по монтажу и эксплуатации**

BINDER GmbH

▶ Address: Post office box 102, 78502 Tuttlingen, Germany ▶ Phone: +49 7462 2005 0
▶ Fax: +49 7462 2005 100 ▶ Internet: http://www.binder-world.com
▶ E-mail: info@binder-world.com ▶ Service Hotline: +49 7462 2005 555
▶ Service Fax: +49 7462 2005 93 555 ▶ Service E-Mail: customerservice@binder-world.com
▶ Service Hotline USA: +1 866 885 9794 or +1 631 224 4340 x3
▶ Service Hotline Asia Pacific: +852 390 705 04 or +852 390 705 03
▶ Service Hotline Russia and CIS: +7 495 988 15 16

Version 07/2022 Art. No. 7001-0269

EU-Konformitätserklärung / EU Declaration of Conformity Déclaration de conformité UE / Declaración de conformidad UE / Декларация соответствия EU

BINDER
Best conditions for your success

EU-Konformitätserklärung / EU Declaration of Conformity / Déclaration de conformité UE / Declaración de conformidad UE / Dichiarazione di conformità UE / Декларация соответствия EU

Hersteller / Manufacturer / Fabricant / Fabricante / Fabbricante / Производитель	BINDER GmbH
Anschrift / Address / Adresse / Dirección / Indirizzo / Адрес	Im Mittleren Ösch 5, 78532 Tuttlingen, Germany
Produkt / Product / Produit / Producto / Prodotto / Продукт	Wasseraufbereitungssystem Water treatment system Système de traitement des eaux Sistema de tratamiento de agua Sistema di trattamento dell'acqua Система очистки воды
Typenbezeichnung / Type / Type / Tipo / Tipo / Тип	BINDER Pure Aqua Service
Art. No. / Art. no. / Réf. / Art. Nº / Art. n. / № арт.	8012-0759

Die oben beschriebenen Produkte sind konform mit folgenden EU-Richtlinien:
The products described above are in conformity with the following EU Directives:
Les produits décrits ci-dessus sont conformes aux directives UE suivantes:
Los productos descritos arriba cumplen con las siguientes directivas de la UE:
I prodotti sopra descritti sono conformi alle seguenti direttive UE:
Продукты, указанные выше, полностью соответствуют следующим EU руководствам:

- **2014/35/EU**
 Niederspannungsrichtlinie 2014/35/EU / Low voltage directive 2014/35/EU / Directive basse tension 2014/35/UE / Directiva sobre baja tensión 2014/35/UE / Direttiva Bassa tensione 2014/35/UE / Директива по низкому напряжению 2014/35/EU
- **2014/30/EU**
 EMV-Richtlinie 2014/30/EU / EMC Directive 2014/30/EU / Directive CEM 2014/30/UE / Directiva CEM 2014/30/UE / Direttiva EMC 2014/30/UE / Директива ЭМС 2014/30/EU
- **2011/65/EU, (EU) 2015/863**
 RoHS-Richtlinien 2011/65/EU und (EU) 2015/863 / RoHS Directives 2011/65/EU and (EU) 2015/863 / Directives RoHS 2011/65/UE et (UE) 2015/863 / Directivas RoHS 2011/65/UE y (UE) 2015/863 / Direttive RoHS 2011/65/UE et (UE) 2015/863 / Директивы RoHS 2011/65/EU и (EU) 2015/863

Die oben beschriebenen Produkte tragen entsprechend die Kennzeichnung CE.
The products described above, corresponding to this, bear the CE-mark.
Les produits décrits ci-dessus, en correspondance, portent l'indication CE.
Los productos descritos arriba, en conformidad, llevan la indicación CE.
I prodotti sopra descritti, conformi a quanto sopra, portano il marchio CE.
Данные продукты в соответствии с изложенным выше маркированы знаком CE.

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BINDER
Best conditions for your success

UKCA Declaration of Conformity

Name and address of manufacturer	BINDER GmbH Im Mittleren Ösch 5, 78532 Tuttlingen, Germany
Name and address of UK Authorised Representative	Comply Express Ltd Unit C2, Coalport House, Stafford Park 1, Telford TF3 3BD
Object of the Declaration	Water treatment system
Type Designation	BINDER Pure Aqua Service
BINDER Art. No.	8012-0759

The Objects of the Declaration described above are in conformity with the relevant UK Regulations and UK Guidelines:

- **Electrical Equipment (Safety) Regulations 2016**
 Statutory Instruments 2016 No. 1101 – Consumer Protection Health and safety
- **Electromagnetic Compatibility Regulations 2016**
 Statutory Instruments 2016 No. 1091 – Electromagnetic Compatibility
- **The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012**
 Statutory Instruments 2012 No. 3032 – Environmental Protection

References of standards and/or technical specifications applied for this Declaration of Conformity, or parts thereof:

S.I. 2016 No. 1101:	EN 61010-1:2010
S.I. 2016 No. 1091:	EN 55011:2016 EN 61000-4-2:2009 EN61000-4-3:2006+A1:2008+A2:2010 EN 61000-4-4:2012 EN 61000-4-5:2014 EN 61000-4-6:2014 EN 61000-4-11:2004
S.I. 2012 No. 3032:	EN IEC 63000:2018

This Declaration is issued under the sole responsibility of the manufacturer.

Tuttlingen	01.06.2022	i.V.	i.V.	
Place	Date	P. Wimmer Vice President	J. Bollaender Director R & D	BINDER GmbH

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Content

1. Safety 29
1.1. Personnel Qualification 29
1.2. Mounting and Operating Instructions 29
1.3. Legal considerations 29
1.3.1. Intellectual property 30
1.4. Structure of the safety instructions 30
1.4.1. Signal word panel 30
1.4.2. Safety alert symbol 31
1.4.3. Pictograms 31
1.4.4. Word message panel structure 31
1.5. General safety instructions on installing and operating BINDER Pure Aqua Service 32
1.6. Material safety data sheet of the cartridge content 33
1.7. Intended use 33

2. Description of the equipment 33

3. Scope of delivery, transportation, storage, and installation 34
3.1. Scope of delivery 34
3.2. Transport and storage 34
3.3. Installation and ambient conditions 35

4. Mounting 36
4.1. Placing the cartridge in the pressure tank 37
4.2. Freshwater connection 38
4.3. Water connection to the climatic chamber 38
4.4. Power connection of the measuring equipment 38
4.5. Overview of the connected water treatment system 39

5. Start up 41
5.1. Rinsing when commissioning a new cartridge 41
5.2. General procedure for commissioning 42

6. Indications on the measuring equipment 43

7. Replacing the cartridge 44

8. Service / contact 44
8.1. Spare parts list 45

9. Performance, lifetime 45
9.1. Water hardness 45
9.2. Maximum flow rate 46

10. Decommissioning 46

11. Disposal 46

11.1. Disposal of the transport packing 46

11.2. Disposal of the cartridge 46

11.3. Disposal of the measuring equipment on the pressure tank cover and the power supply unit in the Federal Republic of Germany 47

11.4. Disposal of the measuring equipment on the pressure tank cover and the power supply unit in the member states of the EU except for the Federal Republic of Germany 47

11.5. Disposal of the measuring equipment on the pressure tank cover and the power supply unit in non-member states of the EU 48

1. Safety

1.1. Personnel Qualification

The device must only be installed, tested, and started up by personnel qualified for assembly, startup, and operation of the device. Qualified personnel are persons whose professional education, knowledge, experience and knowledge of relevant standards allow them to assess, carry out, and identify any potential hazards in the work assigned to them. They must have been trained and instructed, and be authorized, to work on the device.

The device should only be operated by laboratory personnel especially trained for this purpose and familiar with all precautionary measures required for working in a laboratory. Observe the national regulations on minimum age of laboratory personnel.

1.2. Mounting and Operating Instructions

These Mounting and Operating Instructions are part of the scope of delivery. Always keep them handy for reference in the vicinity of the device. If selling the device, hand over the operating manual to the purchaser.

To avoid injuries and damage, observe the safety instructions given with the instructions.

CAUTION

Dangers due to failure to observe the instructions and safety precautions.

Serious injuries and damage to the device.

- Observe the safety instructions in these Mounting and Operating Instructions.
- Follow the operating procedures in these Mounting and Operating Instructions.
- Carefully read the complete Mounting and Operating Instructions of the device prior to installing and using it.
- Keep the Mounting and Operating Instructions for future reference.

Make sure that all persons who use the device and its associated work equipment have read and understood the Mounting and Operating Instructions.

The Mounting and Operating Instructions are supplemented and updated as needed. Always use the most recent version of the Operating Manual. When in doubt, call the BINDER Service Hotline for information on the up-to-dateness and validity of these Mounting and Operating Instructions.

1.3. Legal considerations

These Mounting and Operating Instructions are for informational purposes only. They contain information for correct and safe installing, start-up, operation and decommissioning of the product. Note: the contents and the product described are subject to change without notice.

Understanding and observing the instructions in these Mounting and Operating Instructions are prerequisites for hazard-free use and safety during operation and maintenance. Images are to provide basic understanding. They may deviate from the actual version of the device. The actual scope of delivery can, due to optional or special design, or due to recent technical changes, deviate from the information and illustrations in these instructions. In no event shall BINDER be held liable for any damages, direct or incidental arising out of or related to the use of these instructions.

These Mounting and Operating Instructions cannot cover all conceivable applications. If you would like additional information, or if special problems arise that are not sufficiently addressed in this manual, please ask your dealer or contact us directly, e.g. by phone at the number located on page one of this manual.

Furthermore, we emphasize that the contents of these Mounting and Operating Instructions are not part of an earlier or existing agreement, description, or legal relationship, nor do they modify such a relationship. All obligations on the part of BINDER derive from the respective purchase contract, which also contains the entire and exclusively valid statement of warranty administration and the general terms and conditions, as well as the legal regulations valid at the time the contract is concluded. The statements in this manual neither augment nor restrict the contractual warranty provisions.

1.3.1. Intellectual property

These Mounting and Operating Instructions are protected by copyright. Any unauthorized copying or disclosure to third parties is strictly prohibited. We reserve the right to take legal action and, if necessary, to assert claims for damages in the event of infringement.

Trademark Information: All BINDER trademarks relating to products or service, as well as trade names, logos and product names used on the website, products and documents of BINDER company are trademarks or registered trademarks of BINDER company (including BINDER GmbH, BINDER Inc.) in the U.S. and other countries and communities of states. This includes word marks, position marks, word/figurative marks, design configurations, figurative marks, and design patents.

Patent Information: BINDER products, categories of products, and accessories may be covered by one or more patents and/or utility models in the U.S. and other countries and communities of states. This information is provided to satisfy the virtual patent marking provisions of various jurisdictions, in particular it is intended to serve as notice under 35 U.S.C. § 287(a). Products and services listed on the BINDER website may be sold individually or as part of a combination product. Additional patent applications may also be pending in the U.S. and other countries and communities of states.

Please visit www.binder-world.com for more information.

1.4. Structure of the safety instructions

In these Mounting and Operating Instructions, the following harmonized denominations and symbols indicate dangerous situations following the harmonization of ISO 3864-2 and ANSI Z535.6.

1.4.1. Signal word panel

Depending on the seriousness and probability of the consequences, dangers are identified with a signal word, the corresponding safety color, and if appropriate, the safety alert symbol.

WARNING
Indicates an imminently hazardous situation that, if not avoided, will result in death or serious (irreversible) injury.

CAUTION
Indicates a potentially hazardous situation which, if not avoided, may result in moderate or minor (reversible) injury.

NOTICE
Indicates a potentially hazardous situation which, if not avoided, may result in damage to the product and/or its functions or of a property in its proximity.

1.4.2. Safety alert symbol

Use of the safety alert symbol indicates **risk of injury**.

Observe all measures that are marked with the safety alert symbol in order to avoid death or injury.

1.4.3. Pictograms

Warning signs			
 Electrical hazard	 Danger of frost	 Harmful substances	 Pressurized device
 Explosive atmosphere			
Mandatory action signs		**Prohibition signs**	
 Read instructions	 Environment protection	 Non-potable water	

1.4.4. Word message panel structure

Type / cause of hazard.

Possible consequences.

∅ Instruction how to avoid the hazard: prohibition

➢ Instruction how to avoid the hazard: mandatory action

Observe all other notes and information not necessarily emphasized in the same way, in order to avoid disruptions which could result in direct or indirect injury or property damage.

Please observe the information given on the cartridge label.

1.5. General safety instructions on installing and operating BINDER Pure Aqua Service

IP type of protection: IP 65 (tight against dust and water jet).

NOTICE

Risk of damage by water entering the equipment.

Damage to the equipment.

∅ No liquid must pass by the cable connection into the inside of the measuring equipment.

➢ Ensure dry installation environment.

To avoid any possibility of water damage, ensure that there is a floor drain in the vicinity of where the instrument is located.

Water treated with BINDER Pure Aqua Service is not adapted as drinking water. Drinking the treated water can implicate loss of electrolytes.

⚠ CAUTION

Health hazard by ingestion of non-potable water.

Damages to health.

∅ Do NOT drink processed water.

The cartridge must not freeze since it was filled with water the first time.

NOTICE

Danger of frost.

Damage to the cartridge.

∅ Do NOT expose the cartridge to temperatures below +4 °C / *39.2 °F*.

The connected cartridge is a pressurized device. Do not open. The content may irritate eyes and skin.

⚠ CAUTION

Risk of injury through opening the cartridge by force or damaging it.

Irritation of skin and eyes. Damage to the cartridge and its environment.

∅ Do NOT open the cartridge.

➢ Prevent the cartridge from tipping over.

After using the equipment, close off the water supply at the freshwater supply.

1.6. Material safety data sheet of the cartridge content

	You can order the complete material safety data sheet of the cartridge content at BINDER service.

The following indications refer to the cartridge content. They are only relevant if the content emerges, e.g. due to damage of the cartridge.

Extract from the material safety data sheet of the cartridge content

- Content: Functionalized styrene / divinylbenzene copolymer: Polystyrene sulphonic acid and polyvinyl benzyl dimethyl amine or polystyrene sulphonic acid and polyvinyl benzyl trimethyl ammonium hydroxide.
- Avoid eye and skin contact. Keep away from food and beverages. Eye contact: rinse with plenty of water and seek medical attention. Contact with skin: wash with plenty of water. Ingestion: drink 200-300 ml (half pint) water and seek medical attention
- Do not discharge the cartridge content into drains.

1.7. Intended use

The water treatment system BINDER Pure Aqua Service serves for tap water treatment, preparing water suitable for operating BINDER climatic chambers.

Water treated with the water treatment system BINDER Pure Aqua Service is not suitable for analytical and / or medical purpose. It is no drinking water.

Any use of the device that does not comply with the requirements specified in these Mounting and Operating Instructions shall be considered improper use.

The device is not classified as a medical device as defined by Regulation (EU) No 2017/745.

The device is designed for setting up inside a building (indoor use). The requirements described in the Mounting and Operating Instructions for installation site and ambient conditions (chap. **Fehler! Verweisquelle konnte nicht gefunden werden.**) must be met.

The device does not dispose of any measures of explosion protection.

Contamination of the device by toxic, infectious or radioactive substances must be prevented

2. Description of the equipment

The water treatment system BINDER Pure Aqua Service consists of a pressure tank with measuring equipment to assess water quality, an exchangeable cartridge for water treatment, and the necessary connecting elements.

It serves for tap water treatment.

Dimensions BINDER Pure Aqua Service: height 55 cm / *21.6 in*, width 29 cm / *11.4 in*, depth 26 cm / *10.2 in*

Weight BINDER Pure Aqua Service with transport packing (state of delivery, including cartridge): 20 kg / *44.1 lb*

Weight BINDER Pure Aqua Service, cartridge not filled with water: 18 kg / *39.7 lb*

Weight BINDER Pure Aqua Service, cartridge filled with water: 24 kg / *52.9 lb*

Power consumption of measuring equipment: below 1 W

Technical data of the power supply unit: 100 V-240 V, 50/60 Hz, sec. 12 V, 630 mA

3. Scope of delivery, transportation, storage, and installation

3.1. Scope of delivery

The water treatment system BINDER Pure Aqua Service (Art. no. 8012-0759) consists of:

- Disposable cartridge (Art. no. 6011-0165)
- Pressure tank with handles, cover and measuring equipment with power supply unit 100 V-240 V and water drain for pressure compensation when replacing the cartridge
- Mounting instructions (Art. no. 7001-0269)
- 2 m metal water hose with screwings, for processed water, and 2 gaskets
- Silicon tube for water drain

Figure 1: Scope of delivery BINDER Pure Aqua Service

3.2. Transport and storage

Permissible ambient temperature range during transport and storage:

- Cartridge not filled with water (state of delivery): -10 °C / *14°F* up to +40 °C / *104°F*.
- Cartridge filled with water: +4 °C / *39.2 °F* up to +40 °C / *104°F*.

The cartridge must not freeze since it was filled with water the first time.

	**NOTICE**
	Danger of frost. **Damage to the cartridge.** ∅ Do NOT expose the cartridge to temperatures below +4 °C / *39.2 °F*.

Permissible ambient humidity: 85 % r.h. max., non-condensing

Intermediate storage in a closed and dry room.

3.3. Installation and ambient conditions

Set up BINDER Pure Aqua Service upright on an even surface. The connected cartridge is a pressurized device. The content may irritate eyes and skin.

CAUTION

Risk of injury through opening the cartridge by force or damaging it.

Irritation of skin and eyes. Damage to the cartridge and its environment.

∅ Do NOT open the cartridge.

➢ Prevent the cartridge from tipping over.

Permissible ambient temperature range during operation: +4 °C / *39.2 °F* up to +40 °C / *104°F* (climatic chamber: +18 °C / *64.4°F* to +32 °C / *89.6°F*.)

Permissible ambient humidity: 85 % r.h. max., non-condensing.

To avoid any possible water damage, provide a floor drain at the location of the device. Select a suitable installation site to avoid any consequential damage by splashing water.

To completely separate the unit from the power supply, you must disconnect the power plug. Install the unit in a way that the power supply unit is easily accessible and can be easily disconnected in case of danger.

The unit must not be installed and operated in potentially explosive areas.

WARNING

Explosion hazard through explosive dust or air-solvent mixture in the vicinity of the device

Serious injury or death from burns and / or explosion pressure.

∅ Do NOT operate the unit in potentially explosive areas.

➢ KEEP explosive dust or air-solvent mixtures AWAY from the vicinity of the unit.

4. Mounting

Connections overview


Figure 2: Connections on the pressure tank cover


Figure 3: Connections on the pressure tank cover (mounted)

4.1. Placing the cartridge in the pressure tank

- Remove the pressure tank cover. For this purpose, fix den the pressure tank (place both feet on the black treads on the base of the pressure tank). Press the locking button (inscription “PRESS AND OPEN”) and turn the cover counterclockwise


Figure 4: Removing the pressure tank cover

- Place the cartridge in the pressure tank. Remove the red cover from the cartridge. You can use it for transport of a used cartridge.
- Check the O-ring gasket of the cartridge for correct fit.
- Replace the cover on the pressure tank. For this purpose, fix den the pressure tank (place both feet on the black treads) and put on the cover. The protruding arrow mark must be situated exactly above the arrow of the “INSERT” groove. Press down the cover and turn it clockwise (inscription “CLOSE”),until the locking button locks into place.

Figure 5: Installing the pressure tank cover

4.2. Freshwater connection

Freshwater temperature: not below +5 °C / *41°F* and not above 40 °C / *104°F*.

Supply pressure: 1 bar up to 10 bar

Admissible substances: Max admissible chlorine content: 0.3 mg/l.

- Connect the freshwater hose delivered with the climatic chamber (R3/4“union nut, brass) to the water inlet on the pressure tank cover. The inlet is labeled “IN” and marked with an arrow in the direction of the water flow.
- Connect the other end of the freshwater hose (R3/4” union nut, brass) to a lockable freshwater tap.

4.3. Water connection to the climatic chamber

- Connect the hose for purified water (metal hose with R3/4’’ union nut) to the water outlet on the pressure tank cover (R3/4 in male thread). Insert one of the two gaskets delivered with the water hose. Conduct the hose to the climatic chamber.

Figure 6: Hose for purified water with gaskets at the (state of delivery)

- Connect the hose (metal hose with R3/4’’ union nut) to the freshwater connection at the rear of the climatic chamber. Insert one of the two gaskets delivered with the water hose. Observe the indications on water connection given in the chamber’s operating manual.
- All water connections must be checked for tightness and leaks. This is the operator’s responsibility.

4.4. Power connection of the measuring equipment

- Connect the power supply unit to the measuring equipment
- Prior to connection, check the power supply voltage. A socket with a distribution voltage between 100 V and 240 V / 50/60 Hz is required.

NOTICE

Danger of incorrect power supply voltage due to improper connection.

Damage to the device.

➢ Check the power supply voltage before connection and start-up.

- Only use the original power supply unit according to the above specification.
- IP type of protection: IP 20

	NOTICE
	Risk of damage by water entering the equipment **Damage to the equipment.** ∅ No liquid must pass by the cable connection into the inside of the measuring equipment. ➢ Ensure dry installation environment.

	To completely separate the device from the power supply, you must disconnect the power plug. Install the device in a way that the power plug is easily accessible and can be easily pulled in case of danger.

4.5. Overview of the connected water treatment system


Figure 7: BINDER Pure Aqua Service with connected freshwater hose, hose for purified water and power supply unit

Figure 8: KBF with connected water treatment system BINDER Pure Aqua Service

5. Start up

5.1. Rinsing when commissioning a new cartridge

- Connect the supplied silicon tube to the water drain

Figure 9: Connections on the pressure tank cover, water drain with connected tube

- Introduce the free end of the silicon tube into a receptacle of a volume of at least 10 liters and secure it against moves
- Open the water drain on the pressure tank cover (turn it anti-clockwise)


Figure 10: Closed water drain
(position during operation)


Figure 11: Open water drain
(to replace or rinse the cartridge)

- Slowly turn on the freshwater supply tap and rinse the system with approx. 10 liters of water
- Close the water drain
- All water connections must be checked for tightness and leaks. This is the operator’s responsibility.

5.2. General procedure for commissioning

- Close the water drain on the pressure tank cover
- Slowly turn on the freshwater supply tap

The water treatment system BINDER Pure Aqua Service is now ready for use.

With a new cartridge, the first green lamp on the measuring equipment lights up (chap. 6).

NOTICE

Risk of damage by water entering the equipment.

Damage to the equipment.

∅ No liquid must pass by the cable connection into the inside of the measuring equipment.

➢ Ensure dry installation environment.

Water treated with BINDER Pure Aqua Service is not adapted as drinking water. Drinking the treated water can implicate loss of electrolytes.

CAUTION

Health hazard by ingestion of non-potable water.

Damages to health.

∅ Do NOT drink processed water.

The connected cartridge is a pressurized device. Do not open. The content may irritate eyes and skin.

CAUTION

Risk of injury through opening the cartridge by force or damaging it.

Irritation of skin and eyes. Damage to the cartridge and its environment.

∅ Do NOT open the cartridge.


➢ Prevent the cartridge from tipping over.

After using the equipment, close off the water supply to the measuring equipment at the freshwater tap.

6. Indications on the measuring equipment


Lighting of green lamps

- Very good water quality.


Additional lighting of yellow lamps

- Very good water quality.
- Order a replacement cartridge Art. No. 6011-0165 at BINDER service (chap. 8).

Additional lighting of red lamps

- 1 red lamp: Still adequate water quality. Replace the cartridge.
- 2 or 3 red lamps: Adequate water quality is no longer assured.

Note: With a longer duration of use, the time after which the next lamp lights up diminishes

7. Replacing the cartridge

- Close the freshwater supply tap. Connect the supplied silicon tube to the water drain.

Figure 12: Water drain with connected silicon tube

- Introduce the free end of the tube, which is connected to the water drain into a receptacle of a volume of at least 0.5 liters / *0.018 cu.ft.*
- Open the water drain. Due to the overpressure inside the cartridge, approx. 100 ml water will run off.
- Close the water drain. Open the pressure tank cover and take out the cartridge.
- Close the cartridge for transport with the red cover. You can now dispose of the cartridge (chap. 11.2).
- All water connections must be checked for tightness and leaks. This is the operator's responsibility.

8. Service / contact

Order number for replacement cartridge: **6011-0165**

Contact BINDER Service for orders:

BINDER telephone hotline:	+49 (0) 7462 2005 555
BINDER fax hotline:	+49 (0) 7462 2005 93555
BINDER e-mail hotline:	customerservice@binder-world.com
BINDER service hotline USA:	+1 866 885 9794 or +1 631 224 4340 x3 (toll-free in the USA)
BINDER service hotline Asia Pacific:	+852 390 705 04 or +852 390 705 03
BINDER service hotline Russia and CIS	+7 495 988 15 16
BINDER Internet website	http://www.binder-world.com
BINDER address	BINDER GmbH, post office box 102, 78502 Tuttlingen, Germany

8.1. Spare parts list

Description	Art. No.
Pressure tank cover with measuring equipment and water drain	6014-0029
Power supply unit for measuring equipment 100 V-240 V 50/60Hz sec. 12 V	5020-0047
Handle set, consisting of 2 tank handles and 2 cover handles	6002-0520
Replacement cartridge for BINDER Pure Aqua Service (E2)	6011-0165
EPDM flexible hose, DN 10, 200 cm / *78.7 inches* in length	6008-0250

9. Performance, lifetime

The lifetime of the disposable cartridge depends on the water quality (water hardness) and the water consumption (chap. 9.1).

If the freshwater contains salts due to a preceding water treatment installation or additions of chlorine, the lifetime of the cartridge may be reduced. Max. admissible chlorine content: 0.3 mg/l. Desalinization plants for sea water have no influence.

9.1. Water hardness

Water hardness (freshwater)	Cartridge sufficient for processing the following quantity of water
178.5 mg/l (ppm) = 12.5 English (Clark) degrees = 10.4° grains per U.S. gallon = 1.785 mmol/l	approx. 1300 l
356 mg/l (ppm) = 25 English (Clark) degrees = 20.8° grains per U.S. gallon = 3.57 mmol/l	approx. 350 l
712 mg/l (ppm) = 50 English (Clark) degrees = 41.6° grains per U.S. gallon = 7.14 mmol/l	approx. 390 l

International conversion table for water hardness:

Concentration of alkaline earths ions **mmol/l**	Concentration of calcium carbonate ($CaCO_3$) **mg/l = ppm ***	Concentration of calcium oxide (CaO) **mg/l**	German degrees **°dH**	French degrees **°f**	English degrees / Clark degrees **°E / °Clark**	Concentration of $CaCO_3$ grain per U.S. gallon **gpg ***
1	**100,0**	56.2	5.62	10.0	7.026	5.845
0,017	1.78	**1**	0.10	0.178	0.125	0.104
0,1785	17.85	10	**1**	1.785	1.25	1.04
0,10	10.0	5.6	0.56	**1**	0.700	0.582
0,14	14.2	8.0	0.8	1.428	**1**	0.832
0,17	17.1	9.6	0.96	1.71	1,20	**1**
1.428	**142.8**	**80.0**	**8.0**	**14.28**	**10**	**8.32**
1,785	**178,5**	**100,0**	**10,0**	**17,85**	**12,5**	**10,4**
3,57	**356**	**200,0**	**20,0**	**35,6**	**25**	**20,8**
5,355	**534**	**300,0**	**30,0**	**53,4**	**37,5**	**31,2**
7,14	**712**	**400,0**	**40,0**	**71,2**	**50**	**41,6**

* Both types figure as “American degrees”.

9.2. Maximum flow rate

In order to achieve the specified quality of the treated water, a water flow rate of 5 l / min must not be exceeded.

10. Decommissioning

- Close the water supply at the freshwater tap.
- Open the water drain on the pressure tank cover to release the pressure.
- Disconnect the measuring equipment from the power supply.
- Remove the connections at the freshwater supply and at the climatic chamber.
- Temporal decommissioning: See indications for appropriate storage, chap. 3.2
- Final decommissioning: Dispose of the unit as described in chap. 11.

When taking into operation the system again, please observe the indications given in chap. 5.

11. Disposal

11.1. Disposal of the transport packing

Packing element	Material	Disposal
Transport box with metal clamps	Cardboard	Paper recycling
	Metal	Metal recycling
Top cover	Cardboard	Paper recycling
Edge protection, stuffing	Styropor®	Plastic recycling
	Cardboard	Paper recycling
Bag for operating manual	PE foil	Plastic recycling
Insulating air cushion foil	PE foil	Plastic recycling

If recycling is impossible, all packing parts can also be disposed of with normal waste.

11.2. Disposal of the cartridge

- Dispose the content of the cartridge according to the local and national regulations on a disposal site or in compliance with the synthetic material methods. Dispose to an authorized waste collection point. A controlled combustion is possible. Do not discharge into drains or the environment. In state of delivery (cartridge not filled with water), the material is not classified as hazardous waste (exception: UK).
- EU number for exhausted or saturated material used for the preparation of drinking water or water for industrial use is 19 09 05

 If BINDER Pure Aqua Service was used for treatment of fresh water, you can dispose of the cartridge with normal waste according to applicable local and national regulations (exception: UK).
- EU number for exhausted or saturated material used in waste water treatment plants not otherwise specified is 19 08 06


 If BINDER Pure Aqua Service has not been used for treatment of fresh water and therefore contains harmful impurities, the cartridge must not be thrown in the garbage (household waste). Dispose of the cartridge as hazardous waste according to the applicable local and national regulations.

- In UK, surplus product should be declared a 'Special Waste'. Refer to the 'Control of Pollution (Special Waste) Regulations 1980 - SI 1709'

11.3. Disposal of the measuring equipment on the pressure tank cover and the power supply unit in the Federal Republic of Germany

According to Annex I of Directive 2012/19/EU of the European Parliament and of the Council on waste electrical and electronic equipment (WEEE), BINDER devices are classified as “monitoring and control instruments” (category 9) only intended for professional use“. They must not be disposed of at public collecting points.

The measuring equipment and the power supply unit bear the symbol for the marking of electrical and electronic equipment manufactured / placed on the market in the EU after 13 August 2005 and be disposed of in separate collection according to Directive 2012/19/EU on waste electrical and electronic equipment (WEEE) and German national law for electrical and electronic equipment (Elektro- und Elektronikgerätegesetz, ElektroG). WEEE marking: crossed-out wheeled bin with solid bar under. A significant part of the materials must be recycled in order to protect the environment.

At the end of the device’s service life, have the device disposed of according to the German national law for electrical and electronic equipment (Elektro- und Elektronikgerätegesetz, ElektroG from 20 October 2015, BGBl. I p. 1739) or contact BINDER service who will organize taking back and disposal of the device according to the German national law for electrical and electronic equipment (Elektro- und Elektronikgerätegesetz, ElektroG from 20 October 2015, BGBl. I p. 1739).

NOTICE

Danger of violation against existing law if not disposed of properly.

Failure to comply with applicable law.

∅ Do NOT dispose of BINDER devices at public collecting points.

➢ Have the device disposed of professionally at a recycling company which is certified according to the German national law for electrical and electronic equipment (Elektro- und Elektronikgerätegesetz, ElektroG from 20 October 2015, BGBl. I p. 1739).

or


➢ Instruct BINDER Service to dispose of the device. The general terms of payment and delivery of BINDER GmbH apply, which were valid at the time of purchasing the device.

Certified companies disassemble waste (used) BINDER equipment in primary substances for recycling according to Directive 2012/19/EU. The devices must be free from toxic, infectious or radioactive substances in order to eliminate any health hazards to the employees of the recycling companies.

11.4. Disposal of the measuring equipment on the pressure tank cover and the power supply unit in the member states of the EU except for the Federal Republic of Germany

According to Annex I of Directive 2012/19/EU of the European Parliament and of the Council on waste electrical and electronic equipment (WEEE), BINDER devices are classified as “monitoring and control instruments” (category 9) only intended for professional use“. They must not be disposed of at public collecting points.

The measuring equipment and the power supply unit bear the symbol for the marking of electrical and electronic equipment manufactured / placed on the market in the EC after 13 August 2005 and be disposed of in separate collection according to the Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). WEEE marking: crossed-out wheeled bin with solid bar under.

At the end of the device's service life, notify the distributor who sold you the device, who will take back and dispose of the device according to the Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).

NOTICE

Danger of violation against existing law if not disposed of properly.

Failure to comply with applicable law.

∅ Do NOT dispose of BINDER devices at public collecting points.

➢ Have the device disposed of professionally at a recycling company that is certified according to conversion of the Directive 2012/19/EU into national law.

or

➢ Instruct the distributor who sold you the device to dispose of it. The agreements apply that were agreed with the distributor when purchasing the device (e.g. his general terms of payment and delivery).

➢ If your distributor is not able to take back and dispose of the device, please contact BINDER service.

Certified companies disassemble waste (used) BINDER equipment in primary substances for recycling according to Directive 2012/19/EU. The devices must be free from toxic, infectious or radioactive substances in order to eliminate any health hazards to the employees of the recycling companies.

11.5. Disposal of the measuring equipment on the pressure tank cover and the power supply unit in non-member states of the EU

NOTICE

Danger of violation against existing law if not disposed of properly.

Failure to comply with applicable law. Alteration of the environment.

➢ For final decommissioning and disposal of the device, please contact BINDER service.

➢ Observe the statutory regulations for appropriate, environmentally friendly disposal.

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