

## BINDER Data Logger kits for temperature and humidity

The BINDER Data Logger kits for temperature and humidity can record temperature and humidity data of BINDER equipment. This finely tuned product solution also contains useful accessories for mounting the Logger on the BINDER unit, including cable bushings and a sensor mounting bracket.




### Data Logger kits for temperature and humidity. Technical specifications:

- f Memory capacity: 64000 values
- f Interval of measurement: 1s to 3h
- f Measuring range: -40°C – 70°C / -40°F – 158°F
- f Humidity range: 0 – 100% RH
- f Includes calibration certificate
- f Ambient conditions: -35°C to 55°C / -31°F to 131°F
- f Dimensions (WxHxD): 110x85x35 mm / 4.3x3.4x1.4 inch
- f Power supply: 1 x lithium battery 3.6V, with status display; can be replaced by user

### Data Logger Kit TH 70 | Order No.: 8012-0716




**-40 °C – 70 °C / -40 °F – 158 °F**  **0 – 100 % RH**


#### Scope of supply:

- f Data Logger for temperature and humidity
- f Multi-function sensor (temp/humidity), incl. 2m (6.6 ft) extension cable plus cable bushings and a sensor mounting bracket for the BINDER unit
- f Mounting bracket for attaching the Data Logger onto the BINDER unit

ED/ FD/FED	FP/M	BD/ BF/KB	KBW	MK	MKT	KBF KBF ICH/LQC	KBWF	KMF	MKF
-	-	-	-	-	-	⊕	⊕	-	-

### Data Logger kit TH 70/70 - with additional multi-function sensor | Order No.: 8012-0717



**-40 °C – 70 °C / -40 °F – 158 °F**  **0 – 100 % RH**

#### Scope of supply:

- f Data Logger for temperature and humidity with extra multi-function sensor
- f Multi-function sensor (temp/humidity), incl. 2m (6.6 ft) extension cable plus cable bushings and a sensor mounting bracket for the BINDER unit
- f Mounting bracket for attaching the Data Logger onto the BINDER unit

ED/ FD/FED	FP/M	BD/ BF/KB	KBW	MK	MKT	KBF KBF ICH/LQC	KBWF	KMF	MKF
-	-	-	-	-	-	⊕	⊕	-	-

### Optional: Data Logger software | Order No.: 8012-0821

Software for configuration and evaluation, suitable for all BINDER Data Logger kits, incl. a data cable.



Temperature range sensor



Humidity range sensor

● Commendation

- not available