







Water supply system

| | | Series KBF | Series KBF P | Series KBF LQC | Series KMF | Series KBF-S | |
|-------------------|---|---------------|-----------------|-------------------|---------------|-----------------|---|
| Standard: | Water supply via pressurized, deionized water External | • | • | • | • | – |  |
| Accessory: | BINDER PURE AQUA SERVICE System for preparation or complete desalination of tap water, complete set containing PURE AQUA 300 single-use cartridge, measuring device, and all necessary connecting parts | • | • | • | • | • |  |
| Accessory: | Water supply set Flexible water and waste water management. External, for hanging from the back of the device | • | • | • | • | – |  |
| Accessory: | Water supply set Flexible water and waste water management. External, free-standing | • | • | • | – | – |  |
| Standard: | Convenient water supply Independent water supply via tank | – | – | – | – | • |  |
| Accessory: | Flexible positioning of the water container with magnetic holder, hose adapter for water container with practical quick coupling, and wa- ter stop function. External, with magnetic holder | – | – | – | – | • |  |

• available | – not available