

Water supply system

		Series KBF	Series KBF PRO	Series KBF-S ECO		
Standard:	Water supply via pressurized, deionized water External	•	•	–		
Accessory: Art.No. 8012-0759	BINDER PURE AQUA SERVICE System for preparation or complete desalination of tap water, complete set containing PURE AQUA 300 single-use cartridge, measuring device, and all necessary connecting parts	•	•	•		
Accessory: Art.No. 8012-2592	Water supply set Flexible water and waste water management External, free-standing, with level sensor	•	•	–		
Standard:	Convenient water supply Independent water supply via tank	–	–	•		

• available | – not available

