

## BINDER ALARM KIT

The Alarm Kit comprises the Alarm Kit hardware and the corresponding Alarm Center software



### Function:

Alarm notifications are forwarded to the customer's alarm systems via the Alarm Center software.

The Alarm Center software enables ongoing recordings to be monitored independently of APT-COM™ 4 (without any control options)

- Monitoring of communication between APT-COM™ 4 and unit controller or between APT-COM™ 4 and Alarm Center
- Monitoring of individual tolerance bands for each available parameter

An Ethernet interface for connecting to the Alarm Center monitoring software is provided on the unit.

### Scope of delivery:

- Alarm Kit hardware: Web-IO 4.0 Digital, 2 x in / 2 x out with 1 x switching relay
- The Alarm Center software is included in the APT-COM™ 4 installer
- 24 V/18 W international plug-in power supply unit, with interchangeable European, US and UK adapters on primary side
- 1 m patch cable CAT 5 RJ45 for network connection
- Alarm Kit operating manual

### Technical data

Exterior dimensions (WxHxD): 235 mm x 115 mm x 255 mm  
 Input voltage of plug-in power supply unit: 100 V to 240 V, 50/60 Hz  
 Protection class: IP 20  
 Zero-voltage relay contact: Maximum switching capacity  
 DC: 192 W | 24 V

Article number: 9052-0014

