

Data Logger Kits – temperature and humidity

BINDER Data Logger Kits – make independent recording of temperature and humidity data in the BINDER unit possible. The tailored product solution contains helpful accessories such as the cable access assistance and the sensor fixture, as well as those for mounting the logger to the BINDER unit.



TECHNICAL DATA – DATA LOGGER KITS FOR TEMPERATURE AND HUMIDITY:

- Network data logger
- 8 channels for temperature, humidity, and other parameters
- Storage capacity: 64,000 values
- Measuring interval: 1 s to 3 h
- Measuring range: -50 °C to 100 °C
- Humidity range: 0 to 100% RH
- With calibration certificate
- Ambient conditions: -35 °C to 55 °C
- Dimensions (WxHxD): 104 x 141 x 47 mm
- Power supply: 230 V / 24 V plug-in power supply

Data Logger Kit NET 100 | Order no.: 8012-1837



-50 °C – 100 °C 0 – 100% RH

Scope of delivery:

Data Logger – NET
T/H combination sensor with extension cable (2 m)
Incl. access assistance into the BINDER unit and sensor fixture.
Holder for fixing the Data Logger onto the BINDER unit

ED/FD/ FED	FP/M	BD/ BF/KB	KBW	MK	MKT	KBF/KBF ICH/LQC	KBWF	KMF	MKF
-	-	-	-	-	-	-	-	•	•

Data Logger Kit NET 100/60 – with additional combi sensor | Order no.: 8012-1838



-50 °C – 100 °C 0 – 100% RH

Scope of delivery:

Data Logger – NET with additional combi sensor (T: -40 – 60 °C / RH: 0 – 100%)
T/H combination sensor with extension cable (2 m)
Incl. access assistance into the BINDER unit and sensor fixture.
Holder for fixing the Data Logger onto the BINDER unit

ED/FD/ FED	FP/M	BD/ BF/KB	KBW	MK	MKT	KBF/KBF ICH/LQC	KBWF	KMF	MKF
-	-	-	-	-	-	-	-	•	•

Option: Data Logger Software | Order no.: 8012-0821

Configuration and evaluation software for all BINDER Data Logger Kits, incl. data cable.



Sensor temperature range



Sensor humidity range

• Recommended

- Not available