

Dry-air purge



Individual extras for your simulation chamber

ADVANTAGES FOR CUSTOMERS

- KBF, KMF, and KBWF operational area has been expanded to include operating area C (see image above)
- Can be used in Classic.Line, Avantgarde.Line, KB, KBxx, M, and FP
- Units featuring controlled humidity also control the dry-air purge
- The tried-and-tested Donaldson ULTRAPAC SMART is built in as standard – other manufacturers and models are available on request.



BINDER INDIVIDUAL inquiries
go2binder.com/en-BINDER-INDIVIDUAL