

# Reinforced floor and/or inner chamber



## Individual extras for your simulation chamber

### ADVANTAGES FOR CUSTOMERS

- Depending on the model, each chamber can tolerate loads of up to 300 kg
- The inner chamber reinforcement allows the shelves to bear higher loads – up to 300 kg, depending on the chamber
- Can be used in Classic.Line, KB, KBxx, M, FP, and MKxx

BINDER INDIVIDUAL inquiries  
[go2binder.com/en-BINDER-INDIVIDUAL](https://go2binder.com/en-BINDER-INDIVIDUAL)