

Stacking options for BINDER CO2 incubators. C and CB series

CHAMBER WITH ROLLER BASE	C 150/170			CB (E3/E6)			CB (E7)		
Chamber model	C 150 (E1/E2)	C 170 (E3)	C 170 (E7)	CB 53/60	CB 150/160	CB 210/220	CB 56	CB 170	CB 260
Height of chamber [in]	32.3	34.3	34.3	28.3	36.2	42.1	26.4	34.3	40.2
Height of roller base [in]	4.6	4.6	4.6	4.6	4.6	4.6	4.6	4.6	4.6
Height of chamber on roller base [in]	36.9	38.9	38.9	33.0	40.9	46.8	31.0	38.9	44.8
Art. No. Roller base	9051-0024	9051-0029	9051-0028	9051-0031	9051-0028	9051-0029	9051-0043	9051-0028	9051-0044

Roller base



Roller base with model CB 170



STACK WITH STACKING ADAPTER	C 150/170			CB (E3/E6)			CB (E7)		
Upper chamber model	C 150 (E1/E2)	C 170 (E3)	C 170 (E7)	CB 53/60	CB 150/160	CB 210/220	CB 56	CB 170	CB 260
Lower chamber model	C 150 (E1/E2)	C 170 (E3)	C 170 (E7)	CB 53/60	CB 150/160	CB 210/220	CB 56	CB 170	CB 260
Height of chamber [in]	32.3	34.3	34.3	28.3	36.2	42.1	26.4	34.3	40.2
Height of stacking adapter [in]	3.5	3.5	3.5	3.5	3.5	3.5	3.5	3.5	3.5
Height of two chambers with stacking adapter (in)	68.1	72.0	72.0	60.2	75.9	87.8	56.3	72.0	83.8
Art. No. Stacking adapter	9051-0025	9051-0036	9051-0035	9051-0037	9051-0035	9051-0036	9051-0038	9051-0035	9051-0039


Stacking adapter



Stacking adapter with two C 170



STACK WITH ROLLER BASE AND STACKING ADAPTER	C 150/170			CB (E3/E6)			CB (E7)		
Height of roller base [in]	4.6	4.6	4.6	4.6	4.6	4.6	4.6	4.6	4.6
Height of two chambers with stacking adapter on roller base [in]	72.7	76.7	76.7	64.8	80.6	92.4	60.9	76.7	88.5
Art. No. Roller base	9051-0024	9051-0029	9051-0028	9051-0031	9051-0028	9051-0029	9051-0043	9051-0028	9051-0044



Roller base and Stacking adapter with two C 170



STACK ON STACKING FRAME WITH CASTERS	C 150/170			CB (E3/E6)			CB (E7)		
Upper chamber model	-	-	-	-	CB 150/160	CB 210/220	-	CB 170	CB 260
Lower chamber model	-	-	-	-	CB 150/160	CB 210/220	-	CB 170	CB 260
Height of two chambers on stacking frame [in]	-	-	-	-	80.9	92.7	-	77.0	88.8
Art. No. Stacking frame	-	-	-	-	9051-0020	9051-0023	-	9051-0040	9051-0041

Stack on stacking frame with casters



Stack on stacking frame with casters with model CB 170 und C 170



Further stack combinations can be requested by contacting BINDER INDIVIDUAL: go2binder.com/us-BINDER-INDIVIDUAL