


© Binder GmbH 2017 Diese Zeichnung ist urheberrechtlich geschützt. Alle Rechte vorbehalten.



voltage and frequency changer can be positioned as required

BINDER Im Mittleren Ösch 5
78532 Tuttlingen


MKF 1020 with voltage/frequency
changer

Dokument/
document: **9020-0448** 20

Dateiname/File name: 9020-0448
Druckdatum/printing Date: 22.04.2022

General tolerance ISO 2768-mK	measuring unit inch	Scale 1:25	DIN A3 Sheet 1 of 2
----------------------------------	------------------------	---------------	------------------------

© Binder GmbH 2017 Diese Zeichnung ist urheberrechtlich geschützt. Alle Rechte vorbehalten.



Optional for access port
Ø1,2in, Ø2,0in, Ø3,9in, Ø4,9in leftside

Optional for access port
Ø1,2in, Ø2,0in, Ø3,9in, Ø4,9in rightside

Standard access port Ø2,0in

Standard access port Ø2,0in

Optional for access port
Ø1,2in, Ø2,0in, Ø3,9in, Ø4,9in top

--- inner chamber
 [hatched box] potential position for access ports,
 via BINDER INDIVIDUAL

BINDER Im Mittleren Ösch 5
78532 Tuttlingen

MKF 1020 with voltage/frequency
changer

Dokument/
document: **9120-0448** 20
Dateiname/File name: 9020-0448
 Druckdatum/printing Date: 22.04.2022

	General tolerance	measuring unit	Scale	DIN A3
	ISO 2768-mK	inch	1:20	Sheet 2 of 2